



Xray Database
Version 1



Xray Database
Version 2



Xray Database
Version 3

+fmi_tube_models
+fmi_heat_storage
+fmi_version

-fmi_large_spot_heat_storage
-fmi_small_spot_heat_storage

+tube_model
+last_known_tube_heat
+last_known_tube_heat_epoch

fmi_tube_models	
PK	<u>tube_model</u>
	large_focal_spot_min_power large_focal_spot_max_power small_focal_spot_min_power small_focal_spot_max_power filament_standby_current

fmi_heat_storage	
PK,FK1 PK PK PK	<u>tube_model</u> <u>focal_spot_size</u> <u>exposure_power</u> <u>exposure_timer</u>
	heat_storage

fmi_tubes	
PK	<u>id</u>
FK1	tube_serial_number tube_model last_known_tube_heat last_known_tube_heat_epoch creation_date notes is_active

fmi_generators	
PK	<u>id</u>
	generator_serial_number inverter_serial_number has_z_deflection filament_boost_gain filament_boost_offset creation_date notes is_active

fmi_cooling_curve	
PK,FK1 PK	<u>tube_model</u> <u>timer</u>
	heat_storage

+tube_model

fmi_boost_currents	
PK	<u>id</u>
FK1 FK2	tube_id generator_id tube_voltage tube_current large_focal_spot_boost_current small_focal_spot_boost_current is_modified

-identifier
+generator_serial_number
+inverter_serial_number
+filament_boost_gain
+filament_boost_offset

fmi_boost_currents_baseline	
PK PK	<u>tube_voltage</u> <u>tube_current</u>
	large_focal_spot_boost_current small_focal_spot_boost_current

-boost_current
+large_focal_spot_boost_current
+small_focal_spot_boost_current


fmi_grids	
PK	<u>id</u>
FK1	tube_id <u>spot_size</u> <u>tube_voltage</u> <u>low_grid</u> <u>high_grid</u>
-	
-	
-	
-	

-boost_current
+large_focal_spot_boost_current
+small_focal_spot_boost_current


fmi_grids_baseline	
PK PK	<u>tube_voltage</u> <u>spot_size</u>
	low_grid high_grid

fmi_version	
PK	<u>version</u>
	creation_date notes

Xray Database
Version 4



Xray Database
Version 5



Xray Database
Version 6

